


Galileo TMAtic 5000 **Fully Automated Tissue Arrayer**



www.isenet.it

video: <https://youtu.be/ePd71JhrZPc>

Why Choose Galileo TMAtic

- **HTS Next Generation Fully Automated Tissue Microarrayer Platform.**
- **Proprietary Galileo TMAtic user interface** (based on Galileo TMA CK series platforms)
- **Open architecture** which allows to work with different types of tissue cassettes (Standard, Macro and 96/384 tube microplate in any combination) for molecular analysis.
- **Precise core selection** with automatic Digital 3 points glass slide overlay
- **Cores insertion depth control** (via software)
- **Digital report:** open office (Excel compatible) and XML to interface with commercial Digital Scanners (Aperio, Hamamatsu).
- **Remote SW**, to create the TMA geometry and design in the Pathology Office.
- **Service and post-sales support world-wide.**

**HTS
120-130
cores/hour**

**Needle
cleaning
at each
cycle**

**Needle
integrity
check at
each cycle**

**Exclusive
Technical
Features**

**Tissue
Block
Bar-code
reading**

**Measurement
of Tissue
Blocks
height**

**Control of
extraction of
tissue
cylinder**

Headquarter::

Integrated Systems Engineering srl

Legal Office: Via Volturmo, 80 CEDRI 571; 20816 Brugherio (MB);

Operative Office: **C/O OpenZone**, Via Meucci n. 3 - Edificio B1 - 4 Piano, -

20091 Bresso (MI) - ITALY Tel. **+39.02.36743314**; Cell: +39 3896851496

e.mail: info@isenet.it - PEC: isenet@zucchetticonsulting.it - webpage: www.isenet.it

C.F. 11803980157 - P.IVA. IT03017660964 - Codice SDI: SUBM70(zero)N

Openzone

science oxygen business

Main Office USA:

ISENET-USA LLC

CIC Philadelphia - Cambridge Innovation Center


3675 Market Street, suite 200 - Philadelphia Pa, 19104 - USA

Phone: **+1.215.789.5866**

e.mail: info@isenet-usa.com - web page: www.isenet.it


EIN: 46-1766310 - DUNS: 458743437

cic





Exclusive Features:

- Able to use different types of Tissue Cassettes in any combination: Standard and Macro; 96/384 vials cassettes for further molecular analysis.



- Needle integrity check at each cycle.
- If needle is broken or bent Galileo TMAc will stop, allowing the operator to change the needle and restart the automatic punching cycle.

- Tissue Block height measurement. (5 measurements per block).
- Minimum paraffin height accepted is 2 mm.



- Bar-code reading of table configuration and tissue blocks (displayed on the screen during the construction process).
- Checks if tissue block mounted is part of the design file.


- Automatic punching cycle (120-130 cores/hour) with needle cleaning at every cycle.



- Control sensor if tissue core is inside of the punch needles. This assures that all cores are extracted.
- In case of missing cores, they can be recovered, manually, at the end of the automatic punching cycle.


**FOR MORE INFORMATIONS PLEASE CONTACT YOUR LOCAL ISENET OFFICE
OR SEND AN E-MAIL TO INFO@ISENET.IT
WWW.ISENET.IT**

TMA workflow




TEMPLATE EDIT

Easy preparation of the TMA geometry: definition of TMA block size; selection of needle size and space between cores; stretch function by mouse operation to increase the number of rows and columns; possibility to exclude rows and columns; save templates geometries for future use.



DESIGN

Selection of tissue core position on the TMA block (recipient); possibility to import the Excel File with all Bar-codes and annotation for each donor block as part of the design; manual or automatic position of the cores with possibility of sequential or random positioning; use of different colours to identify different types of tissue.



CONSTRUCTION



Table configuration and tissue block position check (displayed on the screen). Bar-code readings and check if tissue blocks mounted are part of the design file.

Import of High-Definition glass slide image; manual selection of the "area of interest" (with Aperio ScanScope); save into the image folder the "screen shot image" of the marked glass slide.

Precise core selection by choosing:

- **"digital glass slide image overlay"** (select from the "image folder" the marked slide which will be overlaid automatically using the 3P matching function).
- **No Overlay**, select core position directly on tissue block image, If not satisfied, it is possible to use the stretch, Oops or Redo function.

At the end of core selection, start the "Automatic Coring Process".




DIGITAL REPORT (DATA REVIEW)


- Create the final report in Open Office format (Excel compatible) or XML file (to interface with commercial digital scanners)
- The Excel file contains: Bar-code number and image of the tissue cassette (if Bar-codes are not used); annotation (if included); image of the glass slide with the indication of where the cores where extracted.

REMOTE SOFTWARE


Allows multiple user groups to work simultaneously from their PC to define the TEMPLATE and DESIGN, increasing productivity.



Technical Specifications:



- HP Pavilion 24-K 0035NL All-in-One Computer with Windows 10 Operating System or equivalent
- Standard configuration: **35 Std Blocks + 1 cleaning/waste Bin**
- Motorized stage Märzhäuser (SCAN 380x170mm-4 mm) – Speed: 30 mm/sec.
- Controller Märzhäuser TANGO 3 DT 24Vdc & Joystick.
- Cmos Jenoptik Griphax Arktur 2/3 " (8 Mpixel) Camera
- Optics Kowa f 50mm objective autofocus
- Embedded 1D/2D Bar-code reader



- Easy mounting and quick punching needles replacement
- Wide range of punching needles available. (0.6, 1.0, 1.5, 2.0 mm/diam.

TMA cores per Tissue size / punch needles size

Needle size	Standard Cassette (28x42 mm)	Row & Column	Macro Cassette (52x72 mm)	Row & Column
* 0.6 mm/Dia.	640 cores max.	20x32	1.672 cores max.	38x44
* 1.0 mm/Dia.	300 cores max.	15x20	780 cores max.	26x30
* 1.5 mm/Dia.	140 cores max.	10x40	378 cores max.	18x21
* 2.0 mm/Dia.	80 cores max.	8x10	168 cores max.	12x4

CONSUMABLES

PUNCHING NEEDLES SET

0.6



1.0



1.5



2.0 (mm/diam.)



ACCESSORIES

6 Standard Blocks
Holder



3 Standard Blocks
+ 1 Macro Block Holder



6 Standard Blocks
+ 1 Microplate Holder



2 Macro Blocks
Holder



Waste-Cleaning Bin



**FOR MORE INFORMATIONS PLEASE CONTACT YOUR LOCAL ISENET OFFICE
OR SEND AN E-MAIL TO INFO@ISENET.IT
WWW.ISENET.IT**